

INFLUENZA A+B PCR

Collection Instructions

Materials Required

- 1 Sterile nasopharyngeal swab (thinner swab)
- 1 Sterile oropharyngeal swab (thicker swab)
- 1 Sterile virus transport media (collection tube)
- 1 Biohazard bag
- 1 Absorbent sheet
- 1 Parafilm sealing sheet
- 2 Barcode labels


Return supplies:

- UN3373 Biological substance, Category "B" stamped FedEx Clinical pak
- Corrugated return box
- FedEx Return label

Collection Procedures

NASOPHARYNGEAL COLLECTION

1. Wash hands thoroughly.
2. Please have the patient blow their nose to remove excess mucus.
3. Tilt head back 70 degrees, and gently pass the nasal (thinner) swab through the nose until resistance is encountered or the distance is equivalent to that from nostrils to outer opening of the ear.
4. Leave swab in for several seconds to absorb secretions. Slowly remove swab while rotating.
5. Insert the swab into the vial, making sure the swab tip is covered by the liquid in the vial. The tip of the swab must be immersed in liquid.
6. Snap off the swab at the score mark, and secure cap on the vial, making sure no fluid will leak out of the vial.
7. Please keep the vial in an upright position for 10-15 minutes before preparing the sample for shipping back to the lab.



Nasopharyngeal Collection

Safety Information:

- Always wear appropriate protective equipment when handling any clinical specimen. These include lab coats, disposable gloves, and a face mask.
- Do not touch your face when handling specimens. Periodically change gloves and disinfect surfaces with 10% bleach and 70% ethanol.
- Always wash hands thoroughly with soap after handling specimens.